


适配仪器型号	适配仪器图片	0.1ml 8 连管	0.2ml 8 连管	0.1ml 96 孔板	0.2ml 96 孔板	384 孔板	说明
罗氏 LightCycler 480/480 II	 LightCycler480/480II	PCR-0108-NW		PCR-0196-NW (需要辅助器) PCR-0196-HW			
罗氏 LightCycler 96	 LightCycler96	PCR-0108-NW		PCR-0196-NW (需要辅助器) PCR-0196-HW			
罗氏 LightCycler NANO	 LightCycler Nano	PCR-0108-NW PCR-0108-NC					
ABI7300			PCR-0208-NW PCR-0208-NC PCR-0208-WC PCR-0208-WW		PCR-0296-HC PCR-0296-HW PCR-0296-NC PCR-0296-NW		
ABI7500			PCR-0208-NW PCR-0208-NC PCR-0208-WC PCR-0208-WW		PCR-0296-HC PCR-0296-HW PCR-0296-NC PCR-0296-NW		


ABI7500 Fast		PCR-0108-WC PCR-0108-WW PCR-0108-NW PCR-0108-NC		PCR-0196-NC PCR-0196-NW			
ABI Vlla7		PCR-0108-WC PCR-0108-WW PCR-0108-NW (Fast 模块适用)	PCR-0208-WC PCR-0208-WW (普通模块适用)	PCR-0196-NC PCR-0196-NW (Fast 模块适用)	PCR-0296-NC PCR-0296-NW (普通模块适用)		
ABI Stepone plus		PCR-0108-WC PCR-0108-WW PCR-0108-NW		PCR-0196-NC PCR-0196-NW			
ABI Stepone		PCR-0108-WC PCR-0108-WW PCR-0108-NW		48 孔板			
ABI 7900			PCR-0208-WC PCR-0208-WW (配辅助器)		PCR-0296-NC PCR-0296-NW (配辅助器)		

ABI 7900 Fast		PCR-0108-WC PCR-0108-WW PCR-0108-NW (配辅助器)		PCR-0196-NC PCR-0196-NW (配辅助器)		384 孔板 (384 孔模块适用)	
ABI MiniAmp	 MiniAmp		PCR-0208-WC PCR-0208-WW		PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW		
ABI SimpliAmp	 SimpliAmp		PCR-0208-WC PCR-0208-WW		PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW		
ABI Veriti	 Veriti	PCR-0108-WC PCR-0108-WW PCR-0108-NW (Fast 模块适用)	PCR-0208-WC PCR-0208-WW (普通模块适用)	PCR-0196-NC PCR-0196-NW (Fast 模块适用)	PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW (普通模块适用)	384 孔板 (384 孔模块适用)	
ABI ProFlex	 ProFlex		PCR-0208-WC PCR-0208-WW		48 孔板 24 孔板		

ABI QuantStudio 3	 QuantStudio3	PCR-0108-WC PCR-0108-WW PCR-0108-NW (Fast 模块适用)	PCR-0208-WC PCR-0208-WW (普通模块适用)	PCR-0196-NC PCR-0196-NW (Fast 模块适用)	PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW (普通模块适用)	384 孔板 (384 孔模块适用)	
ABI QuantStudio 5	 QuantStudio5	PCR-0108-WC PCR-0108-WW PCR-0108-NW (Fast 模块适用)	PCR-0208-WC PCR-0208-WW (普通模块适用)	PCR-0196-NC PCR-0196-NW (Fast 模块适用)	PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW (普通模块适用)	384 孔板 (384 孔模块适用)	
ABI QuantStudio 6 Flex	 QuantStudio6Flex	PCR-0108-WC PCR-0108-WW PCR-0108-NW (Fast 模块适用)	PCR-0208-WC PCR-0208-WW (普通模块适用)	PCR-0196-NC PCR-0196-NW (Fast 模块适用)	PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW (普通模块适用)	384 孔板 (384 孔模块适用)	
ABI QuantStudio 7 Flex	 QuantStudio7Flex	PCR-0108-WC PCR-0108-WW PCR-0108-NW (Fast 模块适用)	PCR-0208-WC PCR-0208-WW (普通模块适用)	PCR-0196-NC PCR-0196-NW (Fast 模块适用)	PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW (普通模块适用)	384 孔板 (384 孔模块适用)	
ABI QuantStudio 12 Flex	 QuantStudio12KFlex	PCR-0108-WC PCR-0108-WW PCR-0108-NW (Fast 模块适用)	PCR-0208-WC PCR-0208-WW (普通模块适用)	PCR-0196-NC PCR-0196-NW (Fast 模块适用)	PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW (普通模块适用)	384 孔板 (384 孔模块适用)	

伯乐 Bio-rad Opticon 2	 Opticon2	PCR-0108-WC PCR-0108-WW PCR-0108-NW		PCR-0196-NC PCR-0196-NW			
伯乐 Bio-rad CFX Connect	 CFX Connect	PCR-0108-WC PCR-0108-WW PCR-0108-NW (Fast 模块适用)	PCR-0208-WC PCR-0208-WW (普通模块适用)	PCR-0196-NC PCR-0196-NW (Fast 模块适用)	PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW (普通模块适用)		
伯乐 Bio-rad CFX 96/384	 CFX96/384	PCR-0108-WC PCR-0108-WW PCR-0108-NW		PCR-0196-NC PCR-0196-NW		384 孔板 (384 孔模块适用)	
伯乐 Bio-rad CFX 96 Touch	 CFX96 Touch	PCR-0108-WC PCR-0108-WW PCR-0108-NW		PCR-0196-NC PCR-0196-NW			
安捷伦 Agilent/Stratagena MX3000P	 MX3000P		PCR-0208-WC PCR-0208-WW		PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW		


安捷伦 Agilent/Stratagena MX3005P/4000P	 MX3005P/4000P		PCR-0208-WC PCR-0208-WW		PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW		
安捷伦 Agilent/Stratagena MXG8830A	 MXG8830A	PCR-0108-WC PCR-0108-WW PCR-0108-NW		PCR-0196-NC PCR-0196-NW			
耶拿 Analytik Jena qTOWER ³		PCR-0108-WC PCR-0108-WW PCR-0108-NW (0.1ml 模式适用)	PCR-0208-WC PCR-0208-WW (0.2ml 模式适用)	PCR-0196-NC PCR-0196-NW (0.1ml 模式适用)	PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW (0.2ml 模式适用)		推荐 0.1ml 透明系 列
耶拿 Analytik Jena qTOWER3 84/G						适配 384 孔板	
耶拿 Analytik Jena qTOWER ³ auto		PCR-0108-WC PCR-0108-WW PCR-0108-NW (0.1ml 模式适用)	PCR-0208-WC PCR-0208-WW (0.2ml 模式适用)	PCR-0196-NC PCR-0196-NW (0.1ml 模式适用)	PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW (0.2ml 模式适用)		推荐 0.1ml 透明系 列

耶拿 Analytik Jena TOptical		PCR-0108-WC PCR-0108-WW PCR-0108-NW (0.1ml 模式适用)	PCR-0208-WC PCR-0208-WW (0.2ml 模式适用)	PCR-0196-NC PCR-0196-NW (0.1ml 模式适用)	PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW (0.2ml 模式适用)		推荐 0.1ml 透明系 列
Biometra Tadvanced 96 S/SG	 Biometra Tadvanced 96 S/SG		PCR-0208-WC PCR-0208-WW		PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW		首选透 明系列
Biometra Tadvanced 96 /G	 Biometra Tadvanced 96 /G		PCR-0208-WC PCR-0208-WW		PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW		首选透 明系列
Biometra TRIO	 Biometra TRIO		PCR-0208-WC PCR-0208-WW		48 孔板		首选透 明系列
艾本德 Eppendorf Mastercycler X50	 Mastercycler X50	PCR-0108-WC PCR-0108-WW PCR-0108-NW (0.1ml 模式适用)	PCR-0208-WC PCR-0208-WW (0.2ml 模式适用)	PCR-0196-NC PCR-0196-NW (0.1ml 模式适用)	PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW (0.2ml 模式适用)		首推透 明管

艾本德 Eppendorf Mastercycler nexus	 Mastercycler nexus	PCR-0108-WC PCR-0108-WW PCR-0108-NW (0.1ml 模式适用)	PCR-0208-WC PCR-0208-WW (0.2ml 模式适用)	PCR-0196-NC PCR-0196-NW (0.1ml 模式适用)	PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW (0.2ml 模式适用)		首推透 明管
艾本德 Eppendorf Mastercycler nexus X2	 Mastercycler nexus X2	PCR-0108-WC PCR-0108-WW PCR-0108-NW (0.1ml 模式适用)	PCR-0208-WC PCR-0208-WW (0.2ml 模式适用)	48 孔板, 24 孔板	48 孔板 24 孔板		首推透 明管
艾本德 Eppendorf Mastercycler pro/384	 Mastercycler pro/384	PCR-0108-WC PCR-0108-WW PCR-0108-NW (0.1ml 模式适用)	PCR-0208-WC PCR-0208-WW (0.2ml 模式适用)	PCR-0196-NC PCR-0196-NW (0.1ml 模式适用)	PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW (0.2ml 模式适用)	384 孔板 (384 孔模块适用)	首推透 明管
达安 DA 7600	 DA 7600		PCR-0208-WC PCR-0208-WW		PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW		

力康 CG-02/05	 CG-02/05		PCR-0208-WC PCR-0208-WW		PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW		
力康 Trident 960	 Trident960		PCR-0208-WC PCR-0208-WW		PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW		
宏石 SLAN-48P	 SLAN-48P		PCR-0208-WC PCR-0208-WW		48 孔板 24 孔板		
宏石 SLAN-96P	 SLAN-96P		PCR-0208-WC PCR-0208-WW		PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW		

宏石 SLAN-96S	 SLAN-96S		PCR-0208-WC PCR-0208-WW		PCR-0296-NC PCR-0296-NW PCR-0296-HC PCR-0296-HW		
天隆 Gentier 48S/E	 Gentier 48S/E		PCR-0208-WC PCR-0208-WW		48 孔板 24 孔板		
天隆 Gentier 96E	 Gentier 96E		PCR-0208-WC PCR-0208-WW		PCR-0296-HC PCR-0296-HW PCR-0296-NC PCR-0296-NW		
天隆 TL 988	 TL 988	PCR-0108-WC PCR-0108-WW PCR-0108-NW		PCR-0196-NC PCR-0196-NW			

博日 BIOER LineGene Mini	 LineGene Mini		PCR-0208-NW PCR-0208-NC PCR-0208-WC PCR-0208-WW		48 孔板 24 孔板		
博日 BIOER LineGene K Plus	 LineGene K Plus		PCR-0208-NW PCR-0208-NC PCR-0208-WC PCR-0208-WW		PCR-0296-HC PCR-0296-HW PCR-0296-NC PCR-0296-NW		
博日 BIOER LineGene 9600 Plus	 LineGene 9600 Plus		PCR-0208-NW PCR-0208-NC PCR-0208-WC PCR-0208-WW		PCR-0296-HC PCR-0296-HW PCR-0296-NC PCR-0296-NW		
宝生物 TaKaRa TP700	 TP 700		PCR-0208-NW PCR-0208-NC PCR-0208-WC PCR-0208-WW		PCR-0296-HC PCR-0296-HW PCR-0296-NC PCR-0296-NW		

宝生物 TaKaRa TP900	 TP 900		PCR-0208-NW PCR-0208-NC PCR-0208-WC PCR-0208-WW		PCR-0296-HC PCR-0296-HW PCR-0296-NC PCR-0296-NW		
宝生物 TaKaRa TP950	 TP 950		PCR-0208-NW PCR-0208-NC PCR-0208-WC PCR-0208-WW		PCR-0296-HC PCR-0296-HW PCR-0296-NC PCR-0296-NW		